# ORTHOPEDIC FIXATION PINS

- A. Contents**
- B. Material**
- C. Intended use**
  - 1. Indications
  - 2. Contraindications
  - 3. Warnings
  - 4. Adverse reactions
- D. Safety precautions**
- E. Storage, inspection and sterilization Storage**
- F. Cleaning / Disinfection**
- G. Sterilization / Resterilization**

**A. Contents**  
 The package contains one pin.

**B. Material**  
 The pin is made of:  
 • Stainless steel (ASTM F 138 / ASTM F 139 as well as ISO 5832-1).  
 It is supplied non sterile and available in numerous sizes.

**C. Intended use**  
**1. Indications**  
 Fixation pins are intended to perform as a fixation and stabilization unit of bone fractures or as guidance at insertion of pins into the skeletal system.

**2. Contraindications**  
 Do not use the pins in cases of:  
 - States of health, which exclude adequate assistance or healing process, i.e.:  
 • Impairment of blood supply  
 • Inadequate quantity / quality of bone structure  
 • Extreme adiposis  
 • Previous infection  
 • Strong twist or disposition of shank  
 - Mental conditions which exclude participation in rehabilitation programs (Parkinson's, alcohol or drug abuse, etc.)  
 - Major physical activities, adherent with intense percussion, on which the pins are exposed to excessive pressure  
 - Allergic reactions against components  
 The pins have not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. The pins have not been tested for heating or migration in the MR environment. Therefore do not use the pins in a MR environment.

**3. Warnings**  
 - *The device may only be resterilized if it has not been in touch with bodily fluids.*  
 - *Never use or process damaged or defective instruments. Contact your local representative or dealer for repair or replacement. The use of an Instrument for tasks other than those for which they are intended may result in damaged/broken instruments and/or patient injury.*  
 - *Never use pins if the packaging is damaged.*

**4. Adverse reactions**  
 Adverse reactions may include:  
 • Clinical failure (i.e. pain or injury) due to bending, loosening, wear and tear, fracture of pin, loss of fixation, dislocation and/or migration  
 • Pain and/or abnormal sensations due to the presence of the pin  
 • Primary and/or secondary infections  
 • Allergic reactions to pin material  
 • Injury to vessels, nerves and organs  
 • Hematoma and/or impaired wound healing; hemorrhage

**D. Safety precautions**  
 • Prior to use, thoroughly read these Instructions for Use and become familiar with the surgical technique.  
 • Keep the Instructions for Use accessible to all staff.  
 • The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques. Proper surgical performance of the implantation is the responsibility of the operating surgeon.  
 • The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect pin, incorrect operating techniques, the limitations of treatment methods or inadequate asepsis.  
 • Under no circumstances may modular pin components from different suppliers be combined.  
 • Each patient's record shall document the pin used (name, REF number, lot number).  
 • During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed about post- surgical behavioral requirements.  
 • Damage to the weight-bearing structures can give rise to loosening, dislocation and migration, as well as other complications. To ensure the earliest possible detection of such catalysts of pin dysfunction, the pin must be checked periodically post operatively using appropriate techniques.  
 • Never reuse a pin. Although the pin may appear undamaged, previous stresses may have created non-visible damage that could result in pin failure.  
 • Never use pins if the packaging is damaged. A pin with damaged packaging might be damaged and thus may not be used.

**E. Storage, inspection and sterilization storage**  
 The pin is individually packed in protective packaging that is labeled according to its contents.  
 • Always store the pin in the original protective packaging.  
 • Do not remove the pin from the packaging until immediately before use.  
 • Store the pin in normal hospital environmental conditions.

**F. Cleaning / Disinfection**  
 Products delivered in non-sterile condition, must be cleaned, disinfected, and sterilized prior to use. For cleaning and sterilization, remove the pin from its packaging before cleaning, disinfecting and sterilizing the product. A suitable cleaning, disinfection and sterilization process must be applied by the user. Only pH-neutral cleaning agents should be used. The preparation instructions of the respective cleaning and disinfection agent manufacturer must be considered.

**G. Sterilization / Resterilization**  
 Non-sterile pins as well as pins from open packages which have not been implanted must be sterilized by the user. The pins must be sterilized using a process that has been validated by the health care provider - reference ISO 11134, "Sterilization of health care products - Requirements for validation and routine control -Industrial moist heat sterilization" and AAMI TIR 12, "Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: a guide for device manufacturers."

The health care provider must validate that the steam sterilization cycle is effective in their facility since all steam sterilization chambers are unique.

- Pre-vacuum steam sterilization, wrapped items
- 134° C / 273° F, pre-vacuum for 4 minutes and 15 minutes minimum drying time

**WARNING:**  
*Never use or process damaged or defective instruments. Contact your local representative or dealer for repair or replacement. The use of an Instrument for tasks other than those for which they are intended may result in damaged/broken instruments and/or patient injury.*

- The operating surgeon and operating room team must be thoroughly familiar with the operating technique, as well as the range of pins and instruments to be applied. Complete Information on these subjects must be readily available at the workplace.
- The operating procedure must be explained to the patient, and the patient's understanding of the following information must be documented:
  - The patient is aware of the risks associated with general surgery, orthopedic surgery and with general anesthesia.
  - The patient has been informed about the advantages of the pin procedure and about possible alternative treatments.
  - The pin can fail due to excessive load, wear and tear, or infection.
  - The service life of the pin is determined by body weight and physical activity. The pin must not be subjected to overload through extreme strain, or through work-related or athletic activities.
  - Corrective surgery may be necessary if the pin fails.
  - The patient must have their physician carry out follow-up examinations of the pin at regular intervals.

**WARNING:**  
*Never use pins if the packaging is damaged.*

## SYMBOLS

	Caution		Catalog Number	 Integra York PA, Inc. 589 Davies Drive York, PA 17402 USA toll free phone: 866-854-8300 phone: 717-840-9335 fax: 717-840-9347 www.integralife.com/integra-miltex
	Consult Instructions for Use		Lot Number	
	Do Not Reuse		Quantity	
	Manufacturer			

**Rx ONLY** Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

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