Integra® Miltex®

N-Tralig® Intraligamentary Syringe

Rx ONLY CAUTION: Federal (USA) law restricts this device to sale by or on the order of a dentist or physician.

Description

Integra Miltex N-Tralig Intraligamentary Syringe is used to deliver dental anesthetic to the periodontal ligament. The device is supplied non-sterile.

Reference #	Brand	Description
76-55SS	Miltex	N-Tralig Intraligamentary Syringe
76-55MSS	Miltex	N-Tralig Intraligamentary Syringe, Metric
PN2SS	Miltex	N-Tralig Cartridge Barrel (Solid Wall)
PN3SS	Miltex	N-Tralig Cartridge Barrel (Shield Window)
76-56	Miltex	30° Adapter
76-56M	Miltex	30° Adapter, Metric
PN4	Miltex	N-Tralig Barrel Cap
PN4M	Miltex	N-Tralig Barrel Cap, Metric

Indications For Use

Integra Miltex N-Tralig Intraligamentary Syringe is indicated to be used in conjunction with anesthetic needles and cartridges for injection of anesthetic into the periodontal ligament.

Contraindication

- Do not inject into periodontally involved teeth unless performing an extraction.
- Do not use in procedures involving primary teeth with a developing permanent tooth bud.
- Do not use when cementing crowns or when accurate occlusion is a consideration.

Warning

- Do not inject more than 0.2cc of anesthetic per root (one complete squeeze of the handle). Over dosage may result in post-operative pain and inflammation of the tissues.
- Use caution with higher epinephrine concentrations due to higher risk of ischemic necrosis of adjacent tissue and cardiac symptoms.
- Use caution in procedures involving primary teeth.

Adverse Reactions

- Patient may experience pain during administration of the injection.
- Patient may experience tenderness at the injection site after treatment.
- Patient may experience a subjective sensation that the tooth is high in occlusion.

Precautions

- First and last dosage may not deliver the full 0.2 cc of anesthesia. For accurate dosage, discard first dose and continue with second dose.
- The Integra® Miltex® N-Tralig® Intraligamentary Syringe should not be used as an aspirating syringe.
- Infiltration must be carried out over a period of approximately 30 seconds. Faster infiltration might create pressure strong enough to break the carpule.
- A new carpule and needle must be used for each patient.
- Syringe should not be altered or reconfigured for any reason.
- The anesthetic efficacy of the PDL anesthetic injection technique can be unreliable if the needle is not positioned properly.
- To avoid needle sticks, re-cap needle after injection has been given.
- To avoid leakage during an injection, ensure that needle is straight as it penetrates cartridge diaphragm.
- Do not mix metals when sterilizing syringes. Sterilizing syringes with other non-stainless steel instruments will cause damage to syringe and other instruments.
- Use of 30 gauge short or extra short needles is recommended.
- Only use 1.8ml carpule.

Loading Instructions

- Press down on top of the brake disc (A). Grasp dosage rod knob (B) (while keeping brake disc down) and pull dosage rod out as far as possible.
- Unscrew cartridge barrel cap (C).
- Insert anesthetic carpule, rubber or plastic stopper end first, into barrel (D).
- Replace cartridge barrel cap (C).
- Screw 30 gauge short or extra short injection needle onto the cartridge barrel cap (C). This syringe is designed to accommodate disposable needles with either metal or plastic hubs.
- Gently squeeze injection lever (E) until the automatic dosage stop is reached. A partial dosage of anesthetic fluid will be expelled to assure proper functioning.

Instructions For Use

- Administer topical anesthetic to gingival tissue.
- Place the needle against the tooth at an angle towards the tooth and parallel to the long axis of the tooth. Mesial or distal application is preferred, as those areas are better accessed.
- Administer a few drops of anesthesia to further minimize any discomfort. Slowly slide the needle further through the sulcus, approximately 2mm, until it reaches alveolar bone.
- To make sure the needle is correctly placed in the periodontal ligament, a slight pressure should now be applied on the injection lever (E). Firm resistance to the lever should be felt; this back pressure is necessary to ensure entrance to the ligament. If no back pressure is felt, the needle should be further advanced towards the root of the tooth, or introduced into another area of the periodontal ligament.
- Having made certain that the needle is in the periodontal ligament, squeeze injection lever (E) slowly, releasing a complete dose over a period of approximately 30 seconds, which is the proper time required to dispense one dose of 0.2cc. One complete squeeze of the trigger releases one dose sufficient to anesthetize a single rooted tooth. A tooth with multiple roots requires a dosage of 0.2cc per root.
- The infiltrated anesthesia, when properly applied, is effective immediately.

Unloading Instructions

- Press down on top of the break disc (A).
- Pull back the dosage rod knob (B).
- Remove capped needle from the cartridge barrel cap (C).
- Remove cartridge barrel cap (C) from barrel (D)
- Remove the cartridge and discard.

Cleaning And Sterilization Procedures

Before using the instruments clean, dry, visually examine, and sterilize following the procedures below.

Personnel should follow accepted guidelines as recommended in ANSI/AAMI ST79:2006, A1:2008, A2:2009 -*Comprehensive guide to steam sterilization and sterility assurance in health care facilities.*

Cleaning:

Cleaning should occur as soon as possible after instrument is used.

- Rinse device components together including lumens with running, warm tap water (26°C / 79°F to 36°C / 97°F) for a minimum of one (1) minute. Retract and extend the dosage rod (F) during rinsing by pressing down the break disc (A) and moving the dosage rod knob (B) in and out.
- Use soft bristle brush(es) to scrub the outside and inside surfaces of the device components while rinsing the components under running, warm tap water (26°C / 79°F to 36°C / 97°F) for a minimum of one (1) minute. During scrubbing, retract and extend the dosage rod (F) during brushing by pressing down the break disc (A) and moving the dosage rod knob (B) in and out.

Note 1: Never use steel wool, wire brushes, scalpel blades or highly abrasive detergent or cleansers to remove soil as these will damage the instruments protective surface and lead to corrosion.

- Prepare Miltex EZ-Zyme® All-Purpose Enzyme Cleaner solution as per manufacturer’s directions.

Note 2: Other cleaning reagents/detergents are available but have not been validated for use with this syringe.

- Condition (de-gas) ultrasonic cleaner for a minimum of five (5) minutes.
- Separate dissimilar metals prior to cleaning.
- Place device components in the ultrasonic cleaner.
- Sonicate the device components for a minimum of ten (10) minutes.
- Remove the device from the sonicator and rinse device components together including lumens with running, warm tap water (26°C / 79°F to 36°C / 97°F) for a minimum of one (1) minute. Retract and extend the dosage rod (F) during rinsing by pressing down the break disc (A) and moving the dosage rod knob (B) in and out.
- Thoroughly dry the components with a clean, dry, lint-free wipe or cloth. Force air through the lumens of the barrels using filtered air or a syringe until excess water can no longer be visually seen.
- Visually inspect all device components in a well-lit area for cleanliness and dryness.
- Place two sprays of Integra Miltex lubricant under the advancement block (G).
- Wipe off the excess lubricant.

Sterilization

After cleaning & lubricating the reusable instrument, sterilize it using the following procedure.

- Individually single wrap the clean, dry, disassembled instrument in a Self-Seal Sterilization Pouch. Make sure the pouch is sealed.

Notes 3: Gravity sterilization was validated using Fisherbrand CrossTex® Instant Sealing Sterilization Pouches. Pre-Vacuum sterilization was validated using SPS medical Self-Sealing Sterilization Pouches. Other sterilization pouches may be appropriate for the steam method of sterilization but have not been validated for use with this syringe.

- Place pouch in a sterilization chamber using the following parameters to achieve Sterility Assurance Level (SAL) of 10⁻⁴.

Recommended steam sterilization parameter to achieve Sterility Assurance level of 10⁻⁴:

Sterilizer Type	Configuration	Temperature	Exposure Time
Gravity	Wrapped	132°C (270°F)	15 minutes
Pre-Vacuum	Wrapped	132°C (270°F)	4 minutes
Pre-Vacuum	Wrapped	134°C (273°F)	3 minutes

Notes 4: Other sterilization pouches may be appropriate for the steam method of sterilization but have not been validated for use with this syringe.

Storage Conditions

Store suitably packaged and sterilized instruments in a dry, clean, and dust-free environment.

Maintenance Procedures

Improper, ineffective and insufficient maintenance can reduce the life of an instrument and will invalidate the instrument’s warranty.

Protect Instruments: The use of careful preliminary cleaning, use of neutralized pH solutions, adherence to manufacturer’s instructions and visual inspection, will help to keep instruments performing accurately and free of stains.

Certain compounds are highly corrosive to stainless steel and will cause serious damage. Instruments should never be exposed to:

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The following substances should be avoided whenever possible; rinse with copious amounts of water immediately if instruments are inadvertently exposed to any of the following substances:

- Aluminum chloride
- Mercury chloride
- Barium chloride
- Potassium permanganate
- Bichloride of mercury
- Potassium thiocyanate
- Calcium chloride
- Saline
- Carbolic acid
- Sodium hypochlorite
- Chlorinated lime
- Stannous chloride
- Dakin’s solution

Any kind of corrosion will lead to rust on steel. Rust particles can be transferred from one instrument to another, therefore, remove corroding instruments from service to prevent formation of rust on other instruments.

Protect sharp cutting edges and fine working ends of inserts during all maintenance procedures. Avoid loading heavy items on top of delicate and hollow instruments.

Diagnosing Spots and Stains: It is common for instruments to become stained or spotted. Adhering to proper technique during cleaning and sterilizing procedures will prevent most staining occurrences. The following identifies some of the various instrument-related problems practices may encounter.

- Brown Stains: Detergents containing polyphosphates may dissolve copper elements in the sterilizer resulting in brown stains. A dull blue or brown stain is the result of oxidation on the surface.
- Black Stains: Black stains may be the result of contact with ammonia.
- Light or Dark Spots: Spots are often the result of the mineral content in the water used for rinsing, use of non-neutral instrument or an unclean sterilizer chamber.
- Rust Deposits: It is very unlikely for surgical grade steel to rust. Rust colored spots usually appear in localities where water has high iron content.

Returned Goods Policy

Products must be returned in unopened packages with manufacturer’s seals intact to be accepted for replacement or credit unless returned due to a complaint of product defect. Determination of a product defect will be made by Integra. Products will not be accepted for replacement if they have been in the possession of the customer for more than 90 days.

Repairs And Maintenance

Should your instruments require repair or maintenance, contact Integra for return authorization and address. Instruments returned to Integra for repair must have a statement testifying that each instrument has been thoroughly cleaned and sterilized. Failure to supply evidence of cleaning and disinfection will result in a cleaning charge and delayed processing of your instrument repair.

Warranties And Guarantees

In order to ensure warranties and guarantees, instruments in need of repair should be sent to Integra.

Failure to follow these procedures will invalidate instrument’s warranty and can cause instrument to fail.

Product Information Disclosure

INTEGRA AND ITS SUBSIDIARIES (“INTEGRA”) AND MANUFACTURER EXCLUDE ALL WARRANTIES, EXCEPT INTEGRA’S APPLICABLE STANDARD WARRANTY WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. NEITHER INTEGRA NOR MANUFACTURER SHALL BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THIS PRODUCT. NEITHER INTEGRA NOR MANUFACTURER ASSUME NOR AUTHORIZE ANY PERSON TO ASSUME FOR THEM ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS.

Integra, the Integra logo, Miltex and N-Tralig are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. © 2014 Integra LifeSciences Corporation. All Rights Reserved. SYMTRAD.FU Rev. A 10/09/14


Integra York PA, Inc. • 589 Davies Drive, York, PA 17402 USA
866-854-8300 USA • 717-840-2763 outside USA • 717-840-9347 fax
Integrallife.com/integra-miltex

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