

Vantage[®] by Integra[®] Miltex[®]

Nail Nippers

Limit uncertainty with
a simplified cleaning and
processing experience



INTEGRA[®]
LIMIT UNCERTAINTY

Vantage Stainless Steel Nail Nippers

Vantage stainless steel nail and tissue nippers are a comparable alternative to economy grade chrome-plated nippers. Vantage nippers are available in a variety of sizes and jaw configurations and are a perfect solution for limited use and remote patient treatment plans.

Comparison of Features and Benefits

Vantage® Stainless Steel Nippers

- ✓ No need to separate instruments equals a simplified cleaning process
- ✓ Satin finish helps to reduce glare and may provide a better grip during use
- ✓ Double-spring design helps to provide smooth, controlled jaw movement

Chrome-plated Nippers

- ✗ Must be separated and processed separately from dissimilar metals
- ✗ Chrome finish reflects light creating glare during use
- ✗ Single-spring design may provide a more staggered, less controlled jaw movement



versus



To view our complete offering of premium instruments for the Podiatry market visit integralife.com/integra-miltex.

Vantage® by Integra® Miltex®

Nail Nippers



V940200



V940207



V940205



V940215



V940210



V940212

Reference	Description	Jaws	Length
V940200	Nail Nipper, Single Spring	concave	4 5/8" (118 mm)
V940207	Nail Nipper, Barrel Spring	concave	4 5/8" (118 mm)
V940205	Nail Nipper, Barrel Spring	concave	5 1/2" (139 mm)
V940215	Nail Nipper, Double Spring	concave, angled	5 1/2" (140 mm)
V940210	Nail Nipper, Double Spring	concave	5 5/8" (142 mm)
V940212	Nail Nipper, Double Spring	straight	5 5/8" (142 mm)

Due to the hand craftsmanship of the instruments, tolerances in product specifications are subject to minor variations in dimensions and finishes.

Vantage® by Integra® Miltex®

Nail Nippers



V940225 (pictured)
V940226



V940227



V940220



V940245



V940250



V940251

Reference	Description	Jaws	Length
V940225	Nail Nipper, Double Spring	straight	4 5/8" (118 mm)
V940226	Nail Nipper, Double Spring	straight	5" (126 mm)
V940227	Nail Nipper, Double Spring	straight	6" (155 mm)
V940220	Tissue Nipper, Ring Handle	straight	4" (102 mm)
V940245	Tissue Nipper, Double Spring	convex	4" (100 mm)
V940250	Tissue Nipper, Single Spring	convex	4 1/2" (114 mm)
V940251	Tissue Nipper, Double Spring	convex	5" (126 mm)

Due to the hand craftsmanship of the instruments, tolerances in product specifications are subject to minor variations in dimensions and finishes.

For more information or to place an order, please contact:
 Integra ■ 589 Davies Drive, York, PA 17402
 866-854-8300 USA ■ 717-840-2763 outside USA ■ 866-854-8400 fax
integralife.com/integra-miltex

