

MILTEX HI-LIGHT

Z-TYPE HYSTERECTOMY & PARAMETRIUM SCISSORS & FORCEPS

IDEAL FOR VAGINAL AND ABDOMINAL
HYSTERECTOMY PROCEDURES

Tungsten Carbide Inserts

Seamlessly bonded to upper and lower blades to provide strong yet precise cutting sharpness and lasting durability.

Offset Handles

Optimizes viewing area.



Special Jaw Design

Longitudinal grooves and tip serrations ensure a firm, yet atraumatic grip on tissue.

Longer Instrument Length Options

Ideal for deep pelvic, abdominal surgery and obese patients.

The Miltex® Z-Type Hysterectomy & Parametrium Scissors and Forceps are manufactured to traditionally high Miltex craftsmanship standards and made from the highest quality German stainless steel forgings. Miltex® Z-Type Scissors have a precise form and bend to cut tissue lateral to the uterus and are configured to precisely match the Miltex® Z-Type Forceps, which provide a strong, yet atraumatic grip of uterine tissue.

Z-TYPE SCISSORS

- Available in 3 lengths for greater range of procedural use: 9" (22.9 cm), 10.5" (26.7 cm) and 12" (30.5 cm)
- Available in 3 blade configurations: Slightly Curved, Curved, and Fully Curved

Z-TYPE FORCEPS

- Available in 4 lengths for greater range of procedural use: 8.25" (21 cm), 9.5" (24.1 cm), 12" (30.5 cm) and 14" (35.6 cm)
- Available in 4 jaw configurations: Straight, Slightly Curved, Curved, and Angled

See reverse side for details and part numbers.

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.

Miltex
Redefining Excellence

www.miltex.com
toll free 866 854-8300
phone 717 840-9335

1007/.3/B/2

SURG-ZTYPE-HI



+H834SURGZTYPEHI02

Z-Type Scissors

REF	Description
30-2101TC	Z-Type Hysterectomy (Parametrium) Scissors, 9" (22.9 cm), Slight Curve, Carb-N-Sert® Blades
30-2102TC	Z-Type Hysterectomy (Parametrium) Scissors, 9" (22.9 cm), Curved, Carb-N-Sert® Blades
30-2103TC	Z-Type Hysterectomy (Parametrium) Scissors, 9" (22.9 cm), Full Curve, Carb-N-Sert® Blades
30-2105TC	Z-Type Hysterectomy (Parametrium) Scissors, 10-1/2" (26.7 cm), Slight Curve, Carb-N-Sert® Blades
30-2106TC	Z-Type Hysterectomy (Parametrium) Scissors, 10-1/2" (26.7 cm), Curved, Carb-N-Sert® Blades
30-2107TC	Z-Type Hysterectomy (Parametrium) Scissors, 10-1/2" (26.7 cm), Full Curve, Carb-N-Sert® Blades
30-2109TC	Z-Type Hysterectomy (Parametrium) Scissors, 12" (30.5 cm), Slight Curve, Carb-N-Sert® Blades
30-2110TC	Z-Type Hysterectomy (Parametrium) Scissors, 12" (30.5 cm), Curved, Carb-N-Sert® Blades
30-2111TC	Z-Type Hysterectomy (Parametrium) Scissors, 12" (30.5 cm), Full Curve, Carb-N-Sert® Blades

Carb-N-Sert® = Tungsten Carbide



Z-Type Forceps

REF	Description
30-1900	Z-Type Hysterectomy (Parametrium) Forceps, 8-1/4" (21 cm), Straight, Flared Shanks
30-1901	Z-Type Hysterectomy (Parametrium) Forceps, 8-1/4" (21 cm), Slight Curve, Flared Shanks
30-1902	Z-Type Hysterectomy (Parametrium) Forceps, 8-1/4" (21 cm), Curved, Flared Shanks
30-1903	Z-Type Hysterectomy (Parametrium) Forceps, 8-1/4" (21 cm), Angled, Flared Shanks
30-1904	Z-Type Hysterectomy (Parametrium) Forceps, 9-1/2" (24.1 cm), Straight, Flared Shanks
30-1905	Z-Type Hysterectomy (Parametrium) Forceps, 9-1/2" (24.1 cm), Slight Curve, Flared Shanks
30-1906	Z-Type Hysterectomy (Parametrium) Forceps, 9-1/2" (24.1 cm), Curved, Flared Shanks
30-1907	Z-Type Hysterectomy (Parametrium) Forceps, 9-1/2" (24.1 cm), Angled, Flared Shanks
30-1908	Z-Type Hysterectomy (Parametrium) Forceps, 12" (30.5 cm), Straight, Flared Shanks
30-1909	Z-Type Hysterectomy (Parametrium) Forceps, 12" (30.5 cm), Slight Curve, Flared Shanks
30-1910	Z-Type Hysterectomy (Parametrium) Forceps, 12" (30.5 cm), Curved, Flared Shanks
30-1911	Z-Type Hysterectomy (Parametrium) Forceps, 12" (30.5 cm), Angled, Flared Shanks
30-1912	Z-Type Hysterectomy (Parametrium) Forceps, 14" (35.6 cm), Straight, Flared Shanks
30-1913	Z-Type Hysterectomy (Parametrium) Forceps, 14" (35.6 cm), Slight Curve, Flared Shanks
30-1914	Z-Type Hysterectomy (Parametrium) Forceps, 14" (35.6 cm), Curved, Flared Shanks
30-1915	Z-Type Hysterectomy (Parametrium) Forceps, 14" (35.6 cm), Angled, Flared Shanks

